



TollFree (800)423-3527  
Phone (847)885-7805  
Fax (847)885-2596

1107 Tower Road  
Schaumburg, IL 60173  
sales@orthoarch.com

**REF:- STATIC IMMERSION TEST ION RELEASE AVERAGE VALUES FOR BUCCAL TUBES AND BRACKETS**

To whom it may concern,

The following stated ranges of each element and sum for each category can be found below:

<b>Horizon Stainless Steel Bracket System:</b>	( $\mu\text{g}/20$ brackets/7 days)
Fe	< 1,500
Cr	< 200
Ni	< 10
Sum	< 1,710

<b>Fluted Direct Bond Buccal Tubes:</b>	( $\mu\text{g}/8$ tubes/7 days)
Fe	< 2,000
Cr	< 1,000
Ni	< 7,500
Sum	< 10,500

<b>Weldable Buccal Tubes Bonded to Bands:</b>	( $\mu\text{g}/8$ weldable buccal tubes/7 days)
Fe	< 31,000
Cr	< 4,500
Ni	< 2,000
Sum	< 37,500

Safety Data Sheets for the products listed above will be included.

Best Regards,

Doug Bukaty

President/Partner