



ORTHO ARCH

NEW PRODUCT SUPPLEMENT  
2023-2024



SMILES AWAIT



# Aligner Essentials




## Aligner Cases



Our ultra slim aligner cases feature a magnetic closure. Three attractive colors (black, silver, white) with removable silicone inserts to easily clean.

Assorted pack of 12 includes 4 of each

700-650



## Aligner Removers



Contoured hook style end inserts under the aligner for easy removal. Individually packaged.

12 multi-color removes per pack.

700-640



## Chewies

Aligner chewies help seat aligners to work properly and feel more comfortable to the patient. Available in mint scented, bubblegum scented, or unscented.

10 individual packs of 2 in one master bag.



700-660

Mint Scented  
Green (10pairs)



700-661

Bubblegum Scented  
Pink (10pairs)



700-662

Unscented  
White (10pairs)



# Separating Abrasive Strips



Separating Abrasive Handle

#700-509



10 per pack

Separating abrasive strips are used to reduce or contour the shape of the interproximal areas of the teeth. Each surgical grade stainless steel strip is electronically bonded with a durable diamond abrasive to resist stretching, breaking, or corrosion. Strips are .003" thick. Very flexible and tear resistant.

SIZE	FINE	MEDIUM	COURSE	DOUBLE SIDED MEDIUM
4 mm	700-510	700-511	700-512	700-513
6 mm	700-514	700-515	700-516	700-517
8 mm		700-518	700-519	700-520

# Rodium Coated Dental Mirror

Rhodium coating provides the sharpest, brightest images to the smallest detail. #5 24MM Cone socket – Front surface design.



4 per pack

#800-052



#800-051

## Beck Scissors - Curved



#800-170

Light Scissor Curved

Light model cuts aligner materials, bands, crowns, and fine wire. Heavy model also cuts aligner materials, bands, crowns, and heavier wire.



#800-171

Heavy Scissor Curved

## A.J Wilcock Australian



A.J. Wilcock stainless steel wire offers superior hardness and tensile strength for a multitude of applications where exceptional resiliency is desired. .022 special + grade, 25 ft.

#229-408






## Flexible Silicone Photo Contrastor

Enhance your intra oral photographs by blackening out the background. Using contrastors makes the image more presentable to patients, lectures, and publications. Flexible silicone coating is comfortable for the patient and autoclavable.

Photo Contrastor Kit (Set of 3)

#801-040

Anterior




#801-041

Occlusal




#801-042

Buccal Lingual




#801-043

## Photo Glass Mirrors




PALATAL CHILD

#801-050




PALATAL CHILD-TEEN

#801-051




PALATAL ADULT

#801-052



BUCCAL

#801-053



LINGUAL

#801-054

These photographic mirrors offer the finest possible color and image transfer. They are manufactured of 1/8" float glass for distortion-free surface and double-side coated with chromium for maximum reflective qualities. Compatible with autoclave sterilization. Do not place into an ultrasonic unit.

# MOVEMINTS

## Clear Aligner Mints



Movemints helps you enjoy the journey and find relief from the traditional challenges that impact your progress. With the mint that fits, you'll love fresh breath, better hydration, and more comfort – all adding up to delicious progress.

100 pk

#700-619



## New Super Chain

Color	Chain ooo (short)	Chain ooo (long)	Chain ooo (continuous)
Silver	#300-001	#300-002	#300-003
Clear	#300-004	#300-005	#300-006
Pearl	#300-011	#300-021	#300-031
Black	#300-012	#300-022	#300-032
Royal Blue	#300-013	#300-023	#300-033
Ruby Red	#300-014	#300-024	#300-034
Kelly Green	#300-015	#300-025	#300-035
Grape	#300-016	#300-026	#300-036
Bubblegum	#300-017	#300-027	#300-037
Aqua Blue	#300-018	#300-028	#300-038





# Auxiliaries & Attachments



**Crimpable Surgical Slide on Hook - .021 X .025 (10 pack)**

The surgical slide on ball hook crimps to any rectangular wire up to .021" x .025" the ultimate hook crimping plier 800-083N is designed to place and crimp securely.

**#270-218**



**Direct Bond Coil Eyelet (10 pack)**

Bondable closed coil with contoured base.

**#270-217**

**Crimpable Curved Multi Hook (10 pack)**

**#430-019 Right**

**#430-020 Left**



**Bondable Double T Hook (10 pack)**

**#270-219**



double hook

## Patient Kits

Plastic zip pouch with essential oral care products. Disposable Toothbrush, Travel Toothbrush, Dental Floss, Sand Timer, Disposable Mirror, Interdental Brush, Floss Threaders, Patient Relief Wax.



(10 per pack)

**Multi color Assortment  
Purple, Green, Red, Orange**

**#700-209**





ORTHO ARCH

QUALITY ORTHODONTIC  
SUPPLIES & INSTRUMENTS



ORTHO ARCH

Scan the QR code to get started  
with OrthoArch today!

