



ORTHO ARCH

NEW PRODUCT SUPPLEMENT



SMILES AWAIT

Aligner Seaters



Mint Scented
Green (10 pairs)

#700-660



Bubblegum Scented
Pink (10 pairs)

#700-661



Unscented
White (10 pairs)

#700-662

Aligner seaters help seat aligners to work properly and feel more comfortable to the patient. Available in mint scented, bubblegum scented, or unscented.

10 individual packs of 2 in one master bag.

Aligner Cases




Our ultra slim aligner cases feature a magnetic closure. Three attractive colors (black, silver, white) with removable silicone inserts to easily clean.

Assorted pack of 12 includes 4 of each color.

#700-650

Aligner Removers



Contoured hook style end inserts under the aligner for easy removal. Individually packaged. 12 multi-color removers per pack.

#700-640

Separating Abrasive Strips

Separating Abrasive Handle

#700-509



Separating abrasive strips are used to reduce or contour the shape of the interproximal areas of the teeth. Each surgical grade stainless steel strip is electronically bonded with a durable carborundum abrasive to resist stretching, breaking, or corrosion. Strips are .003" thick. Very flexible and tear resistant. (10 pack)

SIZE	FINE	MEDIUM	COURSE	DOUBLE SIDED MEDIUM
4 mm	700-510	700-511	700-512	700-513
6 mm	700-514	700-515	700-516	700-517
8 mm		700-518	700-519	700-520





ARCH-ALIGNERS

ORTHO ARCH

Clear Aligners for a perfect smile.
Arch-Aligners for a predictable clinical outcome.





Significant Price Savings

Best Turnaround Times

Contact us today to explore how you can enhance your aligner offerings for your patients.



FASTRACK

CLASS II APPLIANCE

Achieve an ideal Class I platform at the initial stage of treatment.

SIZE	CAT.No	COLOR CODE	LEFT	RIGHT
15 mm	049-006	Yellow	●	● ●
16 mm	049-007	White	○	○ ○
17 mm	049-008	Red	●	● ●
18 mm	049-009	Blue	●	● ●
19 mm	049-010	Green	●	● ●
20 mm	049-011	Pink	●	● ●
21 mm	049-012	Purple	●	● ●
22 mm	049-013	Black	●	● ●
23 mm	049-014	Yellow	●	● ●
24 mm	049-015	White	○	○ ○
25 mm	049-016	Red	●	● ●
26 mm	049-017	Blue	●	● ●
27 mm	049-018	Green	●	● ●
049-001 Assorted kit of 5 sizes: 23-27 mm				

Selection
▼ Ruler


Fastrack



1st Molar
↔ Midline



Ball and socket with built in stops for controlled rotation & up-righting.



Smooth radiused arm.

Low profile elastic hook.

Anatomically contoured bases for the best fit and secure bond.



ELASTIC PROTOCOL

Level 1: ¼" 6 ½ oz
#049-041

Level 2: ¼" 8 ½ oz
#049-042



D2/Fastrack Tweezer

Reverse action tweezer for a secure grip.

#800-029



Super Key Ligature Ties

Super Key ligature ties have excellent elasticity and rebound to deliver a uniform continuous force. Easy finger hole grip. The 12 most popular colors are color matched with our popular Super Chain. (1,000 ties per pack)

Silver #300-001 Chain(short) #300-002 Chain(long) #300-003 Chain(cont) #300-601 Super Key Ligature Ties	Clear #300-004 Chain(short) #300-005 Chain(long) #300-006 Chain(cont) #300-604 Super Key Ligature Ties	Pearl #300-011 Chain(short) #300-021 Chain(long) #300-031 Chain(cont) #300-611 Super Key Ligature Ties	Black #300-012 Chain(short) #300-022 Chain(long) #300-032 Chain(cont) #300-612 Super Key Ligature Ties
Royal Blue #300-013 Chain(short) #300-023 Chain(long) #300-033 Chain(cont) #300-613 Super Key Ligature Ties	Ruby Red #300-014 Chain(short) #300-024 Chain(long) #300-034 Chain(cont) #300-614 Super Key Ligature Ties	Kelly Green #300-015 Chain(short) #300-025 Chain(long) #300-035 Chain(cont) #300-615 Super Key Ligature Ties	Grape #300-016 Chain(short) #300-026 Chain(long) #300-036 Chain(cont) #300-616 Super Key Ligature Ties
Bubble Gum #300-017 Chain(short) #300-027 Chain(long) #300-037 Chain(cont) #300-617 Super Key Ligature Ties	Aqua Blue #300-018 Chain(short) #300-028 Chain(long) #300-038 Chain(cont) #300-618 Super Key Ligature Ties	Neon Green #300-019 Chain(short) #300-029 Chain(long) #300-039 Chain(cont) #300-619 Super Key Ligature Ties	Lavender #300-020 Chain(short) #300-030 Chain(long) #300-040 Chain(cont) #300-620 Super Key Ligature Ties

Providing better stretch and recovery.

Super Chain

This is a proprietary specialty polyurethane material that has superior memory and very high elastic recovery as compared to conventional polyurethane materials (CPM). 15 foot spool.



Super Key/Chain Combo Kit

#300-600

Color Code Chart

#000-005S

NEW!
More Sizes

INTRAORAL ELASTICS

AMBER LATEX INTRAORAL ELASTICS

Bag Color	Purple	Green	Red	Blue	Orange
Diameter (inch / mm)	1/8 / 3.2mm	3/16 / 4.8mm	1/4 / 6.4mm	5/16 / 7.9mm	3/8 / 9.4mm
(A) Light (2½ oz)	#350-102	#350-202	#350-302	#350-402	#350-502
(B) Medium (4½ oz)	#350-104	#350-204	#350-304	#350-404	#350-504
(C) Heavy (6½ oz)	#350-106	#350-206	#350-306	#350-406	#350-506
(E) Extra Heavy (8 oz)	#350-108	#350-208	#350-308	#350-408	

RAINBOW LATEX INTRAORAL ELASTICS

Bag Color	Purple	Green	Red	Blue	Orange
Diameter (inch / mm)	1/8 / 3.2mm	3/16 / 4.8mm	1/4 / 6.4mm	5/16 / 7.9mm	3/8 / 9.4mm
(B) Medium (4½ oz)	#360-104	#360-204	#360-304	#360-404	#360-504
(C) Heavy (6½ oz)	#360-106	#360-206	#360-306	#360-406	#360-506

NON-LATEX INTRAORAL ELASTICS

Bag Color	Purple	Green	Red	Blue	Orange
Diameter (inch / mm)	1/8 / 3.2mm	3/16 / 4.8mm	1/4 / 6.4mm	5/16 / 7.9mm	3/8 / 9.4mm
(A) Light (2½ oz)	#350-122	#350-222	#350-322	#350-422	#350-522
(B) Medium (4½ oz)	#350-124	#350-224	#350-324	#350-424	#350-524
(C) Heavy (6½ oz)	#350-126	#350-226	#350-326	#350-426	#350-526
(E) Extra Heavy (8 oz)	#350-128	#350-228	#350-328	#350-428	-



50 bags
of 100
elastics



Mini-M Bracket Series™



The Mini-M Stainless Steel bracket offers patient comfort and treatment efficiency with its compact, smooth design. It provides precise tooth movement and easy placement for a streamlined bonding process. Available in .022 Roth and MBT Rx's with an 80# foil mesh pad.

022 MINI-M SS ROTH
CASE 5x5 U&L W/HOOKS

#463-906

022 MINI-M SS MBT
CASE 5x5 U&L W/HOOKS

#463-907

Available as
individual brackets
in packs of 5.

Delta Wing Bracket™

Low-friction appliance system



A unique bracket design that incorporates multiple ligation options for a customized treatment plan. One piece Sure-Lok bracket base design for maximum retention.


.022 DELTA WING
BRACKET CASE (5X5 U&L)

#405-800

Available as
individual brackets
in packs of 5.

Brackets

Straight Wire V-Slot Bracket



V-slot mini twin Stainless Steel bracket system available in .022 Roth Rx. Low profile. 80# foil mesh pad.

022 VERTICAL SLOT ROTH
CASE W/ HK 3,4,5


#472-902

022 VERTICAL SLOT ROTH
CASE W/ HK 3 ONLY

#472-903

Available as
individual brackets
in packs of 5.

Mini Edgewise V-Slot



V-slot mini twin Stainless Steel bracket in an edgewise Rx. Low profile. 80# foil mesh pad. Available in .018 & .022 slot.

.022 MINI EDGEWISE V-SLOT
CASE BRACKET 5X5 U&L

#464-900

.018 MINI EDGEWISE V-SLOT
CASE BRACKET 5X5 U&L

#464-910

Available as
individual brackets
in packs of 5.

Accuwise Edgewise Bracket™



A mini twin edgewise Stainless Steel bracket system. Low profile. 80# foil mesh pad. Available in .018 & .022 slot.

.018 ACCUWISE EDGEWISE
CASE 5X5 U&L

#464-911

.022 ACCUWISE EDGEWISE
CASE 5X5 U&L

#464-912

Available as
individual brackets
in packs of 5.

EPIX MOLAR BAND

- Available loose or prewelded with your favorite buccal and lingual attachment.
- Anatomical design - medium temper that easily adapts to the tooth surface.
- Consistent size graduation for an accurate and secure fit.
- Micro-etched interior ensures maximum retention.
- Contoured rolled edges for maximum patient comfort.
- Permanent laser markings for easy identification (sizes 3-29)



STARTER KIT

Band Size	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	Total
(Mid Range) Right							1	1	2	2	2	2	2	2	2	2	2	2	1	1	1							25
(Mid Range) Left							1	1	2	2	2	2	2	2	2	2	2	2	1	1	1							25

INTRODUCTORY KIT

Band Size	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	Total
(Full Range) Right	1	1	1	1	2	2	2	2	2	2	2	3	3	3	3	3	2	2	2	2	2	2	1	1	1	1	1	50
(Full Range) Left	1	1	1	1	2	2	2	2	2	2	2	3	3	3	3	3	2	2	2	2	2	2	1	1	1	1	1	50
(Mid Range) Right								4	4	4	4	4	5	5	4	4	4	4	4									50
(Mid Range) vLeft								4	4	4	4	4	5	5	4	4	4	4	4									50



Leatherette Band Box

A great way to store and identify your band inventory

#900-175B




Epix Molar Band Tray w/ Lid

A great way to store and identify your band inventory

#900-175

FASTPACK

A convenient chairside dispenser designed to organize and securely hold brackets and buccal tubes, enhancing efficiency during bracket placement.




#000-002F

• Compatible for any bracket system.



LP SERIES BUCCAL TUBES

LP Series 1st and 2nd molar direct bond buccal tubes! These low-profile tubes feature a fluted entrance for easier archwire insertion, a ball hook on both first and second molars, and a contoured foil mesh pad for secure bonding. The second molar tube is shorter than the first for an optimal fit. Available in .018 and .022 MBT Rx only, choose from a quad box (40 tubes) or a pack of 5!

	Description	Torque	Tip	Offset	R/L	Color Code	.018	.022
	First Molar Single Non-Convertible	-14°	0°	8°	U6R	Red	#924-530	#924-540
		-14°	0°	8°	U6L	Green	#924-531	#924-541
		-20°	0°	8°	L6R	Blue	#924-532	#924-542
		-20°	0°	8°	L6L	Yellow	#924-533	#924-543
	Second Molar	-14°	0°	8°	U7R	Red	#924-534	#924-544
		-14°	0°	8°	U7L	Green	#924-535	#924-545
		-10°	0°	8°	L7R	Blue	#924-536	#924-546
		-10°	0°	8°	L7L	Yellow	#924-537	#924-547
		1 st Molar Quad Box - (40 pieces - 10 per quadrant)					#924-538	#924-548
		2 nd Molar Quad Box - (40 pieces - 10 per quadrant)					#924-539	#924-549

FLUTED MOLAR QUAD BOX KITS

Introducing our Quad Box fluted molar kits. A great in office organizer of molar buccal tubes. Available in first and second molar kits of 40 fluted direct bond buccal tubes (10 pieces per quadrant). Roth, MBT and 0 degree sets of 10.

Description	.018	.022
FLUTED 1ST MOLAR QUAD BOX (10 EACH = 40) ROTH	#924-980F	#924-982F
FLUTED 1ST MOLAR QUAD BOX (10 EACH = 40) MBT	#924-981F	#924-983F
FLUTED 2ND MOLAR QUAD BOX (10 EACH = 40) ROTH	#924-984F	#924-986F
FLUTED 2ND MOLAR QUAD BOX (10 EACH = 40) MBT	#924-985F	#924-987F
FLUTED FIRST MOLAR QUAD BOX (10 EACH = 40) EDGEWISE	#924-988F	#924-989F

Micro Direct Bond Tube - No Hook



Description	R/L	Torque	.018	.022
MICRO TUBE - NO HOOK (5)	EDGEWISE	0°	#924-150	#924-151
MICRO TUBE - NO HOOK (50)	EDGEWISE	0°	#924-150B	#924-151B
MICRO TUBE - NO HOOK (5)	UR	-10°	#924-152	#924-156
MICRO TUBE - NO HOOK (5)	UL	-10°	#924-153	#924-157
MICRO TUBE - NO HOOK (5)	LR	-25°	#924-154	#924-158
MICRO TUBE - NO HOOK (5)	LL	-25°	#924-155	#924-159

Weldable TBT & DBT Tubes

Maxillary	Description	R/L	.018	.022
	Triple Convertible Occlusal .050 HG -14T BT	U6R	#942-825	#942-925
		U6L	#942-826	#942-926
	Double Convertible Occlusal .050 HG -14T BT	U6R	#942-827	#942-927
		U6L	#942-828	#942-928

(5 pack)

T-Force Bite Corrector



Class II

Class III

- The T-Force is a fixed intermaxillary appliance for Class II and Class III correction.
- Telescopic tubes with nickel titanium springs deliver a gentle constant force.
- Locks securely to the archwire.
- Does not require laboratory work or patient compliance.

#271-650 Small Size (23-32mm)

#271-651 Standard Size (27-36mm)

#271-652 T-Force locking tool



Includes two interchangeable springs, locking tool and instructions

Eruption Appliance

24K gold plated for nickel sensitive patients, 1" (2.5cm) traction chain 3.5mm diameter, 80 gauge foil mesh contoured base.

Eyelet Style

#270-650

(5 pack)



Loop Style

#270-651

(5 pack)



Direct Bond Coil Eyelet



Bondable closed coil with contoured base.

#270-217
(10 pack)


Crimpable Surgical Slide on Hook



The surgical slide on ball hook crimps to any rectangular wire up to .021" x .025" the ultimate hook crimping plier 800-083N is designed to place and crimp securely.

#270-218
(10 pack)

Bondable Double T Hook



#270-219
(10 pack)



Crimpable Curved Multi Hook

#430-019
Right

#430-020
Left

(10 pack)



Bondable Eyelet Hook

The low-profile bondable eyelet hook is expertly designed with an ideal shape to attach one or two Class II or Class III elastics. A micro-etched, contoured bonding base ensures a superior fit and enhanced bond strength. Works in tandem with aligners, D2 and Fastrack appliances.



#270-210
(10 pack)

Weldable Sheath V-Hook

No distal indent

#270-208L
UL/LR



#270-208R
UR/LL



LIGHT CURE BRACKET ADHESIVE



Light Cure Adhesive Kit
- Push Syringe (20g)

#650-101

Light Cure Adhesive Kit
- Centrix Tips (14g)

#650-102

- 4 x 5g Light Cure Adhesive Paste Syringe
- 1x Light Cure Sealant 7ml
- 1x Liquid Etch 7ml

CHEMI CURE BRACKET ADHESIVE



Chemi Cure Adhesive Kit
- Push Syringe (20g)

#650-201

Chemi Cure Adhesive Kit
- Centrix Tips (14g)

#650-202

- 4 x 5g Chemical Cure Adhesive Paste Syringe
- 1x Light Cure Sealant 7ml
- 1x Liquid Etch 7ml

LIGHT CURE BAND CEMENT



Ultra Band Cement
Blue 5g Syringe

#650-305

- Visible light cure band cement.
- Unlimited working time before curing.
- Fluoride releasing to prevent decalcification.
- Light tensile and comprehensive strength.
- Reduces bond failure and washout.
- Clean removal of band or appliance.



Ultra Band Cement
Natural 5g Syringe

#650-306

DUAL CURE BAND CEMENT



SELF MIX DUAL SYRINGE KITS

Self Mix Band Cement
Dual Syringe
Kit Blue 40g

#650-301

4 x 10g Syringe, 25 x Mixing Tips



Self Mix Band Cement
Dual Syringe
Kit Natural 40g

#650-302

4 x 10g Syringe, 25 x Mixing Tips



INDIVIDUAL COMPONENTS

Self Mix Band Cement
Dual Syringe - Blue (10g)

#650-303

Self Mix Band Cement Dual
Syringe - Natural (10g)

#650-304

Self Mix Band Cement
Replacement Tips 25/pk


#650-307

Photo Glass Mirrors

Photographic mirrors provide exceptional color accuracy and image clarity. Made from 1/8" float glass for a distortion-free surface and double-coated with chromium for superior reflection. Autoclave sterilization compatible; not suitable for ultrasonic units.



PALATAL
CHILD

#801-050




PALATAL
ADULT

#801-052



PALATAL
CHILD-TEEN

#801-051



LINGUAL

#801-053



BUCCAL

#801-054

Flexible Silicone Photo Contrastor

Enhance your intra oral photographs by blackening out the background. Using contrastors makes the image more presentable to patients, lectures, and publications. Flexible silicone coating is comfortable for the patient and autoclavable.



Photo Contrastor Kit (Set of 3)

#801-040



#801-041
Anterior




Occlusal
#801-042



#801-043
Buccal Lingual



K-Loop Distalizing Spring

Developed by Dr. Varun Kalra to distalize the molar in a more bodily fashion.



#226-148

(10 pack)



A.J Wilcox Australian

A.J. Wilcock stainless steel wire offers superior hardness and tensile strength for a multitude of applications where exceptional resiliency is desired. .022 special + grade, 25 ft.

#229-408



Microbrush

Disposable applicators bend easily into any angle for the placement of small controlled amounts of material in areas of limited access. No need to touch the tip to bend, push the tip against a flat surface to bend up to 90 degrees. Regular size tip.



#600-770

4 dispenser tubes of 100 each totaling 400. (100 of each color)

Patient Kits

Plastic zip pouch with essential oral care products: Disposable Toothbrush, Travel Toothbrush, Dental Floss, Sand Timer, Disposable Mirror, Interdental Brush, Floss Threaders, Patient Relief Wax.

(10 per pack)



#700-209

Multi color
Assortment
Purple, Green, Red
and Orange



#700-207

Blue



#700-208

Pink

Beck Scissors - Curved

Light model cuts aligner materials, bands, crowns, and fine wire. Heavy model also cuts aligner materials, bands, crowns, and heavier wire.



#800-170

Light Scissor Curved



#800-171

Heavy Scissor Curved

D2/FASTRACK TWEEZER

Firmly holds bite correctors, such as the D2 and Fastrack appliances for precision placement. It is also a great tool to place the Delta Wing brackets.



#800-029

2MM STEP PLIER

Double sided offset beak creates a right or left handed 2mm step up or step down.



#800-280

Rhodium Coated Dental Mirror

Rhodium coating provides the sharpest, brightest images to the smallest detail. #5 24MM Cone socket - Front surface design.



#800-051



Replacement Heads

#800-052
4 per pack

NEW!

MINI MOLD PLUS-ALIGNER BUTTON

LIGHT-CURE YOUR OWN RESIN ATTACHMENTS IN SECONDS!

- Contoured triangular base design improves retention
- Compatible with most aligners default cutout for buttons
- Reduced button profile for improved comfort
- Mushroom shape button makes it easy to place elastics

Molds may be autoclaved and reused.


Aligner Button Tips and Handle (10 Tips)

#700-132



NEW!

Niti Palatal Expander PRO



The new "streamline design" NPE Pro expands the upper arch, rotates, uprights, and distalizes molars using continuous thermal Niti forces. Simply select the desired size as your final result. Easy to place and remove, with no lab fees. Designed for use with a .036" x .072" horizontal lingual sheath.

Description	Sizes
Kit - 1 each: Sizes 32,34,36,38 and 40mm (5 midrange sizes)	#271-400
26mm Thermal Niti Palatal Expander PRO	#271-401
28mm Thermal Niti Palatal Expander PRO	#271-402
30mm Thermal Niti Palatal Expander PRO	#271-403
32mm Thermal Niti Palatal Expander PRO	#271-404
34mm Thermal Niti Palatal Expander PRO	#271-405
36mm Thermal Niti Palatal Expander PRO	#271-406
38mm Thermal Niti Palatal Expander PRO	#271-407
40mm Thermal Niti Palatal Expander PRO	#271-408
42mm Thermal Niti Palatal Expander PRO	#271-409
44mm Thermal Niti Palatal Expander PRO	#271-410

Thermal Niti Palatal BAR



Weldable lingual sheath

#270-207
(10pk)

Maxillary molar Rotation and Stabilization Appliance. A pre-programmed appliance that provides predictable rotation, torquing and expansion. Used with a .036 x .072 horizontal lingual sheath.

Description	Sizes
Kit - 1 each: Sizes 26-44mm (10 sizes)	#271-420
26mm Thermal Niti Palatal Expander Bar	#271-421
28mm Thermal Niti Palatal Expander Bar	#271-422
30mm Thermal Niti Palatal Expander Bar	#271-423
32mm Thermal Niti Palatal Expander Bar	#271-424
34mm Thermal Niti Palatal Expander Bar	#271-425
36mm Thermal Niti Palatal Expander Bar	#271-426
38mm Thermal Niti Palatal Expander Bar	#271-427
40mm Thermal Niti Palatal Expander Bar	#271-428
42mm Thermal Niti Palatal Expander Bar	#271-429
44mm Thermal Niti Palatal Expander Bar	#271-430

Contact Us:

1-800-423-3527

1107 Tower Rd.
Schaumburg, IL 60173

Customercare@orthoarch.com

www.OrthoArch.com



Scan the QR code to get started with OrthoArch today!



ORTHO ARCH