



Any orders placed over \$1,000.00
get **FREE** Ground Shipping
*U.S. only

Summer SPECIALS

Offer ends July 31, 2026



ORTHO ARCH

*Specials not available online

Ask your rep for
our new product
catalog!



FASTRACK CLASS II APPLIANCE



Reliable

Anatomically contoured bases
for the best fit and secure bond.

Precision

Ball and socket with built in stops for controlled rotation
& uprighting.

Comfortable

Smooth radiused arm.

Identification

Individually Color ID dots

Versatile

Low profile elastic hook.

\$85.00/ea

Per pair 15mm - 31mm

\$420.00/Kit

Sizes 23mm - 27mm



SIZE	*SET (Right & Left)	COLOR ID MARKS	*LEFT (Individual)	*RIGHT (Individual)
15 mm	249-006	Yellow	● 249-006L	●● 249-006R
16 mm	249-007	White	○ 249-007L	○○ 249-007R
17 mm	249-008	Red	● 249-008L	●● 249-008R
18 mm	249-009	Blue	● 249-009L	●● 249-009R
19 mm	249-010	Green	● 249-010L	●● 249-010R
20 mm	249-011	Pink	● 249-011L	●● 249-011R
21 mm	249-012	Purple	● 249-012L	●● 249-012R
22 mm	249-013	Black	● 249-013L	●● 249-013R
23 mm	249-014	Yellow	● 249-014L	●● 249-014R
24 mm	249-015	White	○ 249-015L	○○ 249-015R
25 mm	249-016	Red	● 249-016L	●● 249-016R
26 mm	249-017	Blue	● 249-017L	●● 249-017R
27 mm	249-018	Green	● 249-018L	●● 249-018R
28 mm	249-019	Pink	● 249-019L	●● 249-019R
29 mm	249-020	Purple	● 249-020L	●● 249-020R
30 mm	249-021	Black	● 249-021L	●● 249-021R
31 mm	249-022	Gray	● 249-022L	●● 249-022R
	249-001	Assorted kit of 5 sizes (23-27 mm)		

Recommended accessories



Bondable Eyelet Hook
(10 pack)

#270-210

\$24.95/pack



Elast-O-Loop
(10 pack)

#270-209

\$12.50/pack



**Fastrack Class II
Placement Instrument**

Firmly holds Fastrack
appliance for precision
placement

#800-029

\$19.00



SALES2@ORTHOARCH.COM

WWW.ORTHOARCH.COM

1-800-423-3527



Visit us Now!

D2 APPLIANCE CLASS II CORRECTOR

D2 Appliance | Set of 2 | #271-600

NEW! D2 Appliance - Clear | Set of 2 | #271-600C

D2 Level 1 Amber Elastics | 1/4 in. 61/2 oz. | #350-311

D2 Level 2 Amber Elastics | 1/4 in. 81/2 oz. | #350-312

D2 Level 1 Clear LF Elastics | 1/4 in. 61/2 oz. | #350-313

D2 Level 2 Clear LF Elastics | 1/4 in. 81/2 oz. | #350-314

25 Bags of 100 Elastics



Stainless Steel



FREE D2 Elastics
with orders of 5 or
more.

Stainless Steel SPECIAL PRICING

Buy 1-4 **\$139.00/set**

Buy 5-9 **\$129.00/set**

Buy 10-19 **\$119.00/set**

Buy 20+ **\$109.00/set**

Clear

SPECIAL PRICING

Buy 1-4 **\$145.00/set**

Buy 5-9 **\$135.00/set**

Buy 10-19 **\$125.00/set**

Buy 20+ **\$115.00/set**

D2/FASTRACK TWEEZERS



\$19.00

#800-029

Firmly holds bite correctors, such as the D2 and Fastrack appliances for precise placement.



Can effectively treat Class III cases



Bondable Eyelet Hook

(10 pack)

#270-210

\$24.95/pack

Razor PREMIUM CUTTERS

All cutting edges are made from Tungsten Carbide inserts providing superior strength and sharpness.

Distal End Safety Cutter Long Handle

Maximum Cutting Capacity:
Round Wiresup to .020
Rect. Wiresup to .022 x .028
Braided Wiresup to .016
Length: 147mm

#800-213



Distal End Safety Cutter

Maximum Cut Capacity:
Rect. Wiresup to .022x.028
Round Wiresup to .020
Braided Wiresup to .016
Length: 120mm

#800-214



Hard Wire Cutter

Maximum Cut Capacity:
Hard Wiresup to .020
Length: 119mm

#800-223



Ligature Cutter Miniature

Maximum Cut Capacity:
Soft Wiresup to .015
Length: 120mm

#800-217



MATHIEU NEEDLE HOLDER W/SAFETY LATCH



- Fine Tip **#800-120**
Best for SS Ligatures
and Kobys
- Ultra Fine Tip **#800-121**
Best for Elastic Ligatures

\$59.00/each



Inverted latch will not catch gloves

Buy 5 or more

\$49.00/each

Precision Serrations for Non-Slip Grip



BRACKETS

BOOST YOUR INVENTORY!



HORIZON MINI



NOVA S+



NUSTAR GEN 2 SLB

Buy 5 cases at retail

Get 2 FREE!

Buy 10 cases at retail

Get 5 FREE!

Buy 20 cases at retail

Get 15 FREE!

– More brackets, more savings!

Available now in our new Ready-Set Pack



Archwire Mix & Match Specials




- Nitinol
- Stainless Steel
- Multi-Modulus
- Copper Niti
- Beta T (TMA)

Buy 5 packs at retail **get 2 packs FREE**

Buy 10 packs at retail **get 5 packs FREE**

Buy 20 packs at retail **get 12 packs FREE**

FLUTED OR LP SERIES BUCCAL TUBES



Lowest profile!

\$29.95

Per pack of 5



\$190.00/Kit

LP Quad Box

A great organizer of direct bond molar tubes. Available in 1st and 2nd molars, Roth or MBT Rx. Kit of 40 fluted direct bond tubes (10 pieces per quadrant)

ISORETRACT

NEW! DRY MOUTH RETRACTOR



\$120.00
EACH

**BUY 5 ISORETRACT,
RECEIVE 1 FREE**

SKU #900-250

The ISORETRACT® Mouth Retractor System effectively manages saliva while providing full arch exposure for bonding and restorative procedures.





Regular Glide-Ties

1,000 Modules Per Pack

Glide-Ties™

\$19.00 Per Pack/Spool

31 Amazing Colors to choose from!

Glide-Ties 15' Chain Elastic

Cont.

Short

Long

Strength & Precision Ties!



12 Popular Colors

\$199.00/combo kit

choose from continuous, short or long chain

#300-300 Short Combo Kit
#300-300L Long Combo Kit
#300-300C Cnt Combo Kit

Super Key Ligature Ties

Super Key ligature ties have excellent elasticity and rebound to deliver a uniform continuous force. Easy finger hole grip. The 12 most popular colors are color matched with our popular Super Chain. (1,000 ties per pack)



Ligature Ties
26 Ties Per Stick
1,040 Ties Per Pack
12 color choices

\$23.00
Per pack



\$249.00/
combo kit

choose from continuous, short or long chain

Experience better stretch and recovery.

Color	SKU#
Silver	#300-601
Clear	#300-604
Pearl	#300-611
Black	#300-612
Royal Blue	#300-613
Ruby Red	#300-614
Kelly Green	#300-615
Grape	#300-616
Bubblegum	#300-617
Aqua Blue	#300-618
Neon Green	#300-619
Lavender	#300-620

Super Chain

\$23.00 Per Spool



Superior stretch and recovery.

Color	Short	Long	Continuous
Silver	#300-001	#300-002	#300-003
Clear	#300-004	#300-005	#300-006
Pearl	#300-011	#300-021	#300-031
Black	#300-012	#300-022	#300-032
Royal Blue	#300-013	#300-023	#300-033
Ruby Red	#300-014	#300-024	#300-034
Kelly Green	#300-015	#300-025	#300-035
Grape	#300-016	#300-026	#300-036
Bubblegum	#300-017	#300-027	#300-037
Aqua Blue	#300-018	#300-028	#300-038
Neon Green	#300-019	#300-029	#300-039
Lavender	#300-020	#300-030	#300-040



THERMAL NITI PALATAL BAR

NEW!

\$45.00/EA

\$395.00/kit

Kit - 1 each: Sizes 26-44mm (10 sizes)

#271-420

Weldable lingual sheath

#270-207 (10pk)

NITI PALATAL EXPANDER PRO

NEW!

\$89.00/EA

\$345.00 Kit

KIT: 1 ea: Sizes 32,34,36,38 and 40mm (5 midrange sizes)

#271-400

Streamline design expands the upper arch, rotates, uprights, and distalizes molars with continuous thermal Niti forces. Easy to place and remove, no lab fees. Fits .036" x .072" horizontal lingual sheath.

CRIMPABLE SPLIT STOPS

NEW!

- Fits wires up to .021 x .025
- Crimp with bird beak or Weingart plier
- Dead soft for easy crimping

\$24.95/ per pack

#430-011

10/pk One size

MOLAR BAND SPECIALS

Buy 100 bands at retail

Get 25 % OFF

Buy 200 bands at retail

Get 50 % OFF

All Bands



Bicuspid Plain



Bicuspid Prewelded



Halo Molar Band



Bicuspid Band Tray



Halo Molar Band Tray



Epix Molar Band



ALIGNER CASE



Only
\$35.00
#700-650
Assorted black, white
and silver cases
12 packs

SCISSORS



Light
#800-170
Only
\$49.00



Heavy
#800-171

German S. Steel

Our lightweight aligner scissor cuts and trims a smooth pattern, even in the tightest areas.

ALIGNER REMOVER



Only
\$12.95
#700-640
12 multi-color removers per pack

ALIGNER PLIERS

\$69.00/each

Hole Punch
#800-151

Vertical Indent
#800-152

Horizontal Indent
#800-153


Tear Drop
#800-154

Dimple
#800-155

Eraser
#800-156



ALIGNER SEATERS



10 bags of 2
Mint #700-660
Bubblegum #700-661
Unscented #700-662
Only
\$11.95

MIN MOLD PLUS ALIGNER BUTTON



LIGHT-CURE YOUR OWN RESIN ATTACHMENTS IN SECONDS!

\$55.00

#700-132

Aligner Button Tips and Handle (10 Tips)

CERAMIC DIRECT BOND BUTTONS

#270-202
10pk
Only
\$29.95



SMOOTHEEZ ABRASIVE TRIMMERS

\$23.95/
per pack

10 trimmers
per package.



Gray - Fine

#700-631



Green - Medium

#700-632



Red - Medium
Coarse

#700-633



Brown - Coarse

#700-634

MICRO DIRECT BOND BUTTON



\$29.50

Per pack of 20

#270-206

Low profile for patient comfort.
Anatomical base for better retention.


ORTHO-BOND LIGHT CURE BAND CEMENT

\$22.50
each

Ultra Band Cement
Blue 5g Syringe
#650-305



Ultra Band Cement
Natural 5g Syringe
#650-306



ARCH-MARKERS



#800-446

\$19.95
Per pack of 100

**FREE Dispenser with
6 or more packs**

Thick wax tips upon sturdy,
wooden dowels to provide
precise marking, high
visibility and superior
adhesion. Disposable to
prevent cross contamination.

100 per bag.

\$7.95

#800-447



In-stock & Ready to Ship!



Mini Mold®

LIGHT CURE YOUR OWN RESIN ATTACHMENTS IN SECONDS

MINI MOLD - INTRO KIT



\$69.00/
per kit

#700-101

MINI MOLD -STARTER KIT



\$149.00/
per kit

#700-100M+

MINI MOLD - INDIVIDUALS

\$55.00/per pack



Bite Ramp
#700-117 3mm
#700-117L 5mm
8 tips 1 handle



Button
#700-004
8 tips 1 handle



Cuspid Ramp
#700-119S Pado
#700-119 Adult
8 tips 1 handle



Habit Control
#700-111
8 tips 1 handle



Bracket
#700-011
8 tips 1 handle



Wire Bonder
#700-009 Small
#700-001 Large
8 tips 1 handle



Molar Pad
#700-007
8 tips 1 handle



Tube
#700-123
8 tips 1 handle



Mini Mold Handle
#700-002

Molds can be autoclaved.





ARCH-ALIGNERS
ORTHO ARCH

PERFECT SMILE.
PREDICTABLE RESULTS!

ARCH ALIGNERS ARE COMFORTABLE, EFFECTIVE AND DISCRETE.

Treatment Planning

Clear-Arch Aligners are custom-crafted to perfectly fit each patient's smile.

Clear-Arch 10

*Hybrid/Very Simple/
Ortho Relapse Cases*

10 Stages per arch

Clear-Arch 20

Simple Cases

20 Stages per arch

Clear-Arch 30


Moderate Cases

30 Stages per arch

Clear-Arch Unlimited

Complex Cases

Unlimited Stages



MCS - Multi-layer clear aligner material

**BUY 4 CLEAR-ARCH
UNLIMITED CASES
GET ONE FREE
CLEAR-ARCH 10 CASE!**



Net Price can not be combined with the above promotion

Configuration	No. of Aligners per Arch	Cost per Case	Special
Clear Arch 10	10	\$695	\$660
Clear Arch 20	20	\$995	\$945
Clear Arch 30	30	\$1,295	\$1,230
Clear Arch Unlimited	Unlimited	\$1,495	\$1,420
Retainers 2 Sets	2 U/L	\$195	\$185
Retainers 3 Sets	3 U/L	\$295	\$280

Contact us today to explore how you can enhance your aligner offerings for your patients.

1-800-423-3527

service@orthoarch.com

www.OrthoArch.com



Scan Me to Learn More!